## **INSTRUCTIONS FOR HAIR ANALYIS**

- 1. Use pubic hair only in spite of what the lab form says, unless a child or doctor has asked you to use had hear. If you use head hair, follow instruction sheet in lab envelop.
- 2. Use stainless steel scissors to cut hair.
- 3. Cut hair as close to skin as possible.
- 4. Length of hair should not exceed 1 <sup>1</sup>/<sub>2</sub> inches. It may be necessary to trim back hair length first.
- 5. Place all hair cutting s over a clean piece of paper. You will need 2 packed teaspoons (this is typically all your pubic hair).
- 6. Place cuttings into small plastic bag. Make sure all blanks are filled in on bag.
- 7. Place this bag into envelope addressed to lab, fill in their form and mail. Do not send them money, we will bill you for \$70.00 directly.